

Omni-C™ Technology


Delivers SNPs to Genome Conformation In a Single Assay

Powered by a sequence-independent endonuclease, Omni-C™ technology represents the next-evolution in Hi-C assays delivering unbiased long-range reads and highly uniform sequence coverage.

Gain a comprehensive view of the genome from SNPs to indels, SNVs to phasing, and loops to topologically associated domains (TADs), all from a single NGS library.

ACHIEVE SHOTGUN-LIKE COVERAGE

Omni-C libraries contain uniform, shotgun-like coverage enabling downstream applications sensitive to variant detection.



Coverage distribution of different Hi-C libraries (colored lines) compared to a shotgun library (dotted line).

CAPTURE LONG-RANGE INFORMATION

Omni-C libraries offer greater long-range information compared to other Hi-C approaches. The result is greater support for chromatin 3-D interactions such as TADs and loops enabling comprehensive genome conformation mapping.

A.


| Library Type | % cis <1 kbp | % cis >1 kbp | % unique molecules at 300 M read pairs |
|----------------|--------------|--------------|--|
| Omni-C | 5% | 95% | 73% |
| Hi-C Multi-RE | 20% | 80% | 69% |
| Hi-C Single-RE | 34% | 66% | 50% |

Comparison of Omni-C libraries with multi-restriction enzyme (RE) and single-RE Hi-C libraries.

A. Proportion of long-range reads for each library type

B. Contact matrices at 4 kb resolution from 800 M total reads per library. Loops are circled with the number of raw reads supporting each contact indicated.

B.



To place an order or for more information:

visit us at www.dovetailgenomics.com or send an email to info@dovetail-genomics.com

GENETIC VARIANTS AND PHASED GENOMES FROM A SINGLE LIBRARY

The shotgun-like characteristics of Omni-C data enable genotyping and haplotype phasing.

| Library Type | True Positive | False Positive | False Negative | Sensitivity | Precision |
|---------------|------------------|----------------|----------------|--------------|--------------|
| Shotgun | 2,696,291 | 9,270 | 6,814 | 99.7% | 99.7% |
| Omni-C | 2,666,339 | 20,081 | 36,766 | 98.6% | 99.3% |
| RE Based Hi-C | 2,387,235 | 33,554 | 315,870 | 88.3% | 98.6% |


Homozygous SNP calls sampled at 300 M read pairs.

| Library Type | #Variants Phased | % Heterozygous SNPs phased | Switch Error Rate | Largest Phase Block Size | #Chromosomes phased end to end |
|---------------|------------------|----------------------------|-------------------|--------------------------|--------------------------------|
| Shotgun | 2,229,492 | 81.44% | 0.0036 | 28.2 kb | 0 |
| Omni-C | 2,299,248 | 84.00% | 0.0100 | 248.01 Mb | 23 |
| RE Based Hi-C | 1,986,467 | 72.60% | 0.0357 | 247.95 Mb | 23 |

GM12878 haplotype assembly metrics using high confidence heterozygous SNPs from the Illumina Platinum Genome.

EXPLORE THE LANDSCAPE OF STRUCTURAL VARIANTS

Omni-C libraries capture large structural (>1 Mbp) variant information such as indels, inversions and chromosomal translocations.




Validated large SVs found in breast cancer cell line HCC1187 representing an inversion, a deletion and inter-chromosomal translocation

EXTRACT MORE DATA FROM EVERY SEQUENCING RUN

Omni-C technology captures genome-wide topology at single nucleotide resolution in a single library prep.

Therefore, in addition to being ideal for the study of 3-D genome conformation, it is suitable for applications traditionally addressed by whole genome sequencing.



To place an order or for more information:

visit us at www.dovetailgenomics.com or send an email to info@dovetail-genomics.com